

## FEATURES

- $\geq 95$  PFE
- Helps protect against odor and dust
- 3-layer protection
  - Outer Layer: Non-woven fabric
  - Middle Layer: Melt-blown fabric
  - Inner Layer: Soft for your skin and absorbs the mist you breathe out
- Eco-friendly and adjustable plastic clip for a proper fit
- Disposable and maintenance-free
- Color: blue
- Main Material: non-woven fabric + melt-blown fabric + non-woven fabric
- Latex Free

## Face Mask – Single Use



Model #: MSK-F1CO


## APPROVALS AND STANDARDS

- FDA, product code QKR
- Health Canada, Class I, PPE

## MATERIAL COMPOSITION

- Material: polypropylene

### 3 Layers of Protection



## PRODUCT USE INFORMATION

### HOW TO USE:

- 1) Remove mask from bag. Unfold mask and apply to face from bridge of nose down to jaw.
- 2) Lightly press along nose piece to fit.
- 3) Hang elastic straps over both ears.
- 4) Adjust to fit mask to cover from bridge of nose to jaw.
- 5) Dispose after use or if it becomes soiled. Single use only. Do not reuse.

## CAUTION

- Caution: This product is not intended for use except as indicated. Please refer to [www.HoMedics.com/learn/face-masks](http://www.HoMedics.com/learn/face-masks) for additional information.



## MASK PRODUCT SPECIFICATIONS

Product Dims	6.89" x 3.74"
Product Weight	.0006 lbs
Shelf Life	5 years from manufacture

## MSK-F1CO-10B PRODUCT SPECIFICATIONS

Package Dims	7.13" x 3.90" x 0.87"
Case Dims.	8.11" x 5.53" x 7.75" EST
Cube Dims.	0.20 ft <sup>3</sup>
Package Weight	0.19 lbs. EST
Case Weight	2.53 lbs. EST
Units per Pouch	10
Boxes in Master	12
UPC	0 31262 09972 3

## MSK-F1CO-20B PRODUCT SPECIFICATIONS

Package Dims	7.13" x 3.90" x 1.38"
Case Dims.	8.58" x 7.68" x 9.53" EST
Cube Dims.	0.31 ft <sup>3</sup>
Package Weight	0.26 lbs EST
Case Weight	4.02 lbs. EST
Units per Box	20
Boxes in Master	12
UPC	0 31262 09973 0

## MSK-F1CO-50B PRODUCT SPECIFICATIONS

Package Dims	7.13" x 3.90" x 3.11"
Case Dims.	14.88" x 8.58" x 10.39" EST
Cube Dims.	0.69 ft <sup>3</sup>
Package Weight	0.49 lbs. EST
Case Weight	7.49 lbs. EST
Units per Box	50
Boxes in Master	12
UPC	0 31262 09974 7